

**BPM-350 USB digital microscope****Introduction**

BPM-350 USB digital microscope provides powers from 20× and 300× with 5.0MP Image sensor. It is ideally suited for Medical, Industrial Inspection, Engineering, educational and Science Applications to examine coins, stamps, rocks, relics, insects, plants, skin, gems, circuit boards, various materials and many other objects. With the included software, you can observe the magnified images, capture video, take snapshots and do measurement with Windows Win7, Win 8, Win 10 32bit&64 bit, Mac OS X 10.5 or above Operation system.

**Feature**

1. Easy to use portable digital microscope with USB output interface.
2. 20× and 300× Magnification.
3. 5.0 Mega Pixels CMOS sensor.
4. Manual Focus from 10mm to 50mm.
5. LED illumination with 8pcs LED lamps, brightness adjustable.
6. Photo management software, Auto White Balance and Exposure.
7. Powered by USB 2.0 cable, 5V DC.
8. Metal stand to hold the microscope steady.
9. Multi-language support: English, Germany, Spanish, French, Russian.
10. Operating System: Windows Win7, Win 8, Win 10 32bit&64 bit, Mac OS X 10.5 or above.

**Application**

BPM-350 USB digital microscope can be used for: Hobbyists, Educators, Medical Labs, Industrial Inspection, Engineering Applications, Teachers, Students, Science Applications, Doctor's Offices, Police Agencies, Government

Testing, and general use by consumers. It is ideally suited for examining solid objects such as coins, stamps, rocks, relics, insects, plants, skin, gems, circuit boards, various materials, and many other objects.

## Specification

Image Sensor	5.0 Mega Pixels CMOS sensor (true resolution)
Image Capture Resolution	5M, 3M, 1.3M, VGA, QVGA
Video Capture Resolution	5M, 3M, 1.3M, VGA, QVGA
Focus Range	Manual Focus From 10mm to 500m
Frame Rate	Max. 30f/s under 600 Lus Brightness
Magnification Ratio	20x to 300x
Video Format	AVI format
Still Image Format	JPEG or BMP
Light Source	8 LED (adjustable by control wheel)
PC Interface	USB 2.0
Power Supply	5V DC from USB port
Operation System	Windows Win7, Win 8, Win 10 32bit&64 bit, Mac OS X 10.6 or above
OSD Language	English, German, Spanish, French, Russian
Bundled Software	MicroCapture with measurement & calibration function
Size	110mm(L) x33mm(R)
Weight	90g
Package Content	Digital microscope with USB cable, Professional stand, User Manual, CD with Software
Packing Info	Gift box, 12pcs/carton, 11.8kgs/carton, carton size: 45x36x43cm

## Sample Image



Sweater



Currency